



Anti Human HSP90 $\alpha\beta$ Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZHSP20-EX
Clone No.:	Hyb-K3720A The splenic lymphocytes from BALB/c mouse immunized with E.coli ricombinant human Heat-Shock Protein (HSP) 90 β were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-K3720A) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with common region of Human HSP 90α and Human HSP 90β (Amino acid residues 291-304)
Purify:	Ammonium sulfate precipitation
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.4mL [1mg/mL]
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

