



Anti PreS1 Monoclonal Antibody (Mouse)

| | |
|---------------------|--|
| Catalog No.: | SIM-2ZPS41-EX |
| Clone No.: | Hyb-T0606 The splenic lymphocytes from BALB/c mouse immunized with HBV-PreS1 range synthetic peptide were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-T0606) were screened. |
| Source: | BALB/C Mouse implantation ascites |
| Specificity: | This antibody specifically reacts with HBV-PreS1 antigen |
| Purify: | Ammonium sulfate precipitation |
| Subclass: | IgG1 |
| Buffer: | Physiological salt solution containing 0.1% Sodium Azide |
| Quantity: | 0.5mL [0.5mg/mL] |
| Storage: | Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles. |

For research use only; not for use as a diagnostic.

