



For research use only. Not for clinical diagnosis.

Catalog No. 2ZTHBR2

THBR2

Code	2ZTHBR2
Description	Mouse IgM, Purified
Source	Mouse ascites
Application	This product inhibits nonspecific reaction caused by human anti-mouse antibody (HAMA) in immunoassays using mouse monoclonal antibodies. The recommended concentration is 0.01-0.1 mg/mL for ELISA, however, the optimum condition should be determined for each assay.
Purification	Ammonium sulfate precipitation
Subclass	IgM
Form	Saline (containing 0.1% NaN ₃)
Protein Concentration	10mg/mL
Volume	2mL
Storage	Maintain at 2-8 °C (Below -20 °C for long term storage) Avoid repeated freeze/thaw cycles



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619



人と科学のステキな未来へ

コスモ・バイオ株式会社

For research use only. Not for clinical diagnosis.

Catalog No. 2ZTHBR2

THBR2

Code	2ZTHBR2
Description	マウス IgM 精製品
Source	マウス腹水
Application	マウス抗体を用いたイムノアッセイにおける、マウスに対する異好性抗体（HAMA:Human Anti-Mouse Antibody）による非特異反応を抑制します。 サンドイッチ測定系における標準添加濃度：10～100 μ g/mL あらかじめ各測定系での最適条件を検討してください。
Purification	硫酸アンモニウム塩析精製
Subclass	IgM
Form	0.1%アジ化ナトリウムを含む生理的食塩水溶液
Protein Concentration	10mg/mL
Volume	2mL
Storage	2～8℃ 3ヶ月 -20℃ 3年間（凍結融解の繰り返しは避けてください）



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619