



Anti Peptide YY (PYY)(Human) Serum

Cat. No. YII-Y072-EX Lot No. 37140301

Description: This antiserum was raised in a rabbit by immunization with a carrier free synthetic PYY (human) peptide. The product vial contains 50 μ L of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry, or any other immunoreaction with PYY.

Immunogen: Synthetic PYY (human), carrier free **Host:** Rabbit

Amino Acid Sequence of PYY (human)¹⁾:

YPIKPEAPGE DASPEELNRY YASLRHYLNL VTRQRY-NH₂

Product Form: Lyophilized unpurified serum **Size:** 50 μ L

Reconstitution: Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN₃ 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

Suggested Working Dilution Range: 1:3,000-10,000 (final dilution ~1:70,000) for radioimmunoassay;

1: 1,000- 9,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

Specificity (based on radioimmunoassay): PYY (human) 100%, PYY (19-36)(human) <0.001%, PYY (1-19) (human) <0.001%, PYY (rat) 80%, PYY (13-36)(rat) <0.1%, PP (human) 0%, NPY (human) <0.001%

Positive Control (immunohistochemistry): Rat ileum

Species Tested: Rat

REFERENCES:

1) K. Tatemoto, I. Nakano, et al., Isolation and primary structure of human peptide YY. Biochemistry Biophysics Research Communications 157:713-718, 1988

FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



2480-1 Awakura, Fujinomiya-shi,
Shizuoka 418-0011, JAPAN
Phone: +81-544-22-2771
Fax: +81-544-22-2770
URL : www.yanaihara.co.jp
E-mail : ask@yanaihara.co.jp