



**Anti Calcitonin Gene-Related Peptide (α-CGRP) (Mouse, Rat) Serum**  
**Cat. No. YII-Y340-EX      Lot No. 41100224**

**Description:** This antiserum was raised in a rabbit by immunization with a carrier free synthetic α-CGRP (mouse, rat) peptide. The product vial contains 50 μL of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with CGRP peptides.

**Immunogen:** Synthetic α-CGRP (mouse, rat), carrier free **Host:** Rabbit

**Amino Acid Sequence of α-CGRP (rat)<sup>1,2)</sup>:**

SCNTATCVTH RLAGLLSRSG GVVKDNFVPT NVGSEAF

**Product Form:** Lyophilized unpurified serum **Size:** 50 μL

**Reconstitution:** Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN<sub>3</sub> 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

**Suggested Working Dilution Range:** 1:2,000~ 6,000 (final dilution ~1:42,000) for radioimmunoassay;

1: 1,000~5,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

**Specificity** (based on radioimmunoassay): α-CGRP (mouse, rat) 100%

**Positive Control** (immunohistochemistry): Rat colon

**Species Tested:** Rat

**REFERENCES:**

- 1) S.G. Amara, V. Jonas et al., Alternative RNA processing in calcitonin gene expression generates mRNAs encoding different polypeptide products. *Nature* 298:240-244,1982
- 2) S.G. Amara, J.L. Arriza et al., Expression in brain of a messenger RNA encoding a novel neuropeptide homologous to calcitonin gene-related peptide. *Science* 229:1094-1098,1985

**FOR RESEARCH LABORATORY USE ONLY**

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

**Distributor**



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMA BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

**Manufacturer**



株式会社  
矢向原研究所

Yanaihara Institute Inc.  
2480-1 Awakura, Fujinomiya-shi,  
Shizuoka 418-0011, JAPAN  
Phone: +81-544-22-2771  
Fax: +81-544-22-2770  
URL : [www.yanaihara.co.jp](http://www.yanaihara.co.jp)  
E-mail : [ask@yanaihara.co.jp](mailto:ask@yanaihara.co.jp)