



Anti Calcitonin Gene-Related Peptide (α -CGRP) (Mouse, Rat) Serum **Cat. No. YII-Y340-EX Lot No. 41100224**

Description: This antiserum was raised in a rabbit by immunization with a carrier free synthetic α -CGRP (mouse, rat) peptide. The product vial contains 50 μ L of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with CGRP peptides.

Immunogen: Synthetic α -CGRP (mouse, rat), carrier free **Host:** Rabbit

Amino Acid Sequence of α -CGRP (rat)^{1,2}:

SCNTATCVTH RLAGLLSRSG GVVKDNFVPT NVGSEAF

Product Form: Lyophilized unpurified serum **Size:** 50 μ L

Reconstitution: Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN₃ 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

Suggested Working Dilution Range: 1:2,000~ 6,000 (final dilution ~1:42,000) for radioimmunoassay;

1: 1,000~5,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

Specificity (based on radioimmunoassay): α -CGRP (mouse, rat) 100%

Positive Control (immunohistochemistry): Rat colon

Species Tested: Rat

REFERENCES:

- 1) S.G. Amara, V. Jonas et al., Alternative RNA processing in calcitonin gene expression generates mRNAs encoding different polypeptide products. Nature 298:240-244,1982
- 2) S.G. Amara, J.L. Arriza et al., Expression in brain of a messenger RNA encoding a novel neuropeptide homologous to calcitonin gene-related peptide. Science 229:1094-1098,1985

FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

