



Anti ACTH (21-39) (Mouse, Rat) Serum

Cat. No. YII-Y351-EX Lot No. 16190424

Description: This antiserum was raised in a rabbit by immunization with a keyhole limpet hemocyanin (KLH) conjugate of synthetic ACTH (21-39) (mouse, rat) peptide. The product vial contains 50 µL of the titled antiserum obtained by lyophilizing its 0.001M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with the C-terminal portion of ACTH (mouse, rat).

Immunogen: Synthetic ACTH (21-39) (mouse, rat)-KLH conjugate **Host:** Rabbit

Amino Acid Sequence of ACTH (21-39) (mouse, rat)¹⁾:
SYSMEHFRWG KPVGKKRRPV KVYPNVAENE SAEAFPLEF

Product Form: Lyophilized unpurified serum **Size:** 50 µL

Reconstitution: Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN₃ 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

Suggested Working Dilution Range: 1:1,000(final dilution ~1:7,000) for radioimmunoassay; 1:1,000-4,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

Specificity (based on radioimmunoassay):
ACTH (21-39) (mouse, rat) 100%, ACTH (mouse, rat) 100%

Positive Control (immunohistochemistry): Rat pituitary gland

Species Tested: Rat

REFERENCES:

1) J. Drouin, M. Chamberland et al., Structure of the rat pro-opiomelanocortin (POMC) gene. FEBS Letters. 193:54-58, 1985

FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



2480-1 Awakura, Fujinomiya-shi,
Shizuoka 418-0011, JAPAN
Phone: +81-544-22-2771
Fax: +81-544-22-2770
URL : www.yanaihara.co.jp
E-mail : ask@yanaihara.co.jp