



Anti Urocortin 2 (Mouse) Serum

Cat. No. YII-Y363-EX Lot No. 1344170121

Description: This antiserum was raised in a rabbit by immunization with porcine thyroglobulin (pTG) conjugate of synthetic urocortin 2 (mouse) peptide. The product vial contains 50mL of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5 mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with urocortin 2 (mouse).

Immunogen: Synthetic Urocortin 2 (mouse)-pTG conjugate

Host:Rabbit

Amino Acid Sequence of Urocortin 2 (mouse)¹⁾:

VILSLDVPIG LLRILLEQAR YKAARNQAAT NAQILAHV-NH₂

Product Form: Lyophilized unpurified serum Size: 50mL

Reconstitution: Reconstitute the product with 0.5 mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN₃ 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing- thawing should be avoided.

Suggested Working Dilution Range: 1:8,000 for enzyme immunoassay (EIA); 1:5,000 for immunohistochemistry . Optimal dilution should be determined by each laboratory for each application.

Specificity(based on EIA): Urocortin 2 (mouse) 100%, urocortin 2 (7-38) (mouse) 80.52%, urocortin 2 (14-38) (mouse) 24.23%, urocortin 2 (21-38) (mouse) 0.21%, urocortin 2 (26-38) (mouse) 0.41%, urocortin 2 (32-38) (mouse) 0.19%, urocortin 2 (1-13) (mouse) 0.12%, urocortin 1 (mouse, rat) 0%, urocortin 1 (human) 0%, urocortin 2 (rat) < 1%, urocortin 3 (mouse) 0%, ACTH (human) 0%, ACTH (mouse, rat) 0.61%, CRF (human, mouse, rat) 0%, CRF (17-41) (human, mouse, rat) 0%

Positive Control (Immunohistochemistry): Hypothalamus

Species Tested(Immunohistochemistry): Mouse

Related Peptides:

Urocortin 1 (mouse, rat) DDPPLSIDLT FHLRLTLEL ARTQSQRERA EQNRIIFDSV-NH₂

Urocortin 3 (mouse) FTLSLDVPTN IMNILFNIDK AKNLRAKAAA NAQLMAQI-NH₂

Related Antiserum:Anti Urocortin 1 (mouse, rat) Serum Y362

Anti Urocortin 3 (mouse) Serum Y364

REFERENCE:

1) TM. Reyes, K. Lewis et al.,Urocortin II: a member of the corticotropin-releasing factor (CRF) neuropeptide family that is selectively bound by type 2 CRF receptors. *Proc. Natl. Acad. Sci. U. S. A.*, **98**: 2843-2848, 2001

FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



2480-1 Awakura, Fujinomiya-shi,
Shizuoka 418-0011, JAPAN
Phone: +81-544-22-2771
Fax: +81-544-22-2770
URL : www.yanaihara.co.jp
E-mail : ask@yanaihara.co.jp