



Anti Humanin Serum

Cat. No. YII-Y440-EX **Lot No. 1282160726**

Description: This antiserum was raised in a rabbit by immunization with a keyhole limpet hemocyanin (KLH) conjugate of synthetic humanin peptide. The product vial contains 50 mL of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with humanin.

Immunogen: Synthetic humanin-KLH conjugate

Host: Rabbit

Amino Acid Sequence of Humanin ^{1) 2)}:

MAPRGFSCLL LLTSEIDLVP KRRA

Product Form: Lyophilized unpurified serum

Size: 50 mL

Reconstitution: Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN₃ 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing- thawing should be avoided.

Suggested Working Dilution Range: 1:8,000 for enzyme immunoassay (EIA); 1:5,000 for immunohistochemistry . Optimal dilution should be determined by each laboratory for each application.

Specificity (based on EIA): Humanin (5-24) 35.0%, humanin (9-24) 1.9%, humanin (14-24) 0.18%, humanin (19-24) 0.18%, humanin (9-18) 0%, rattin-38 0%, rattin-24 0%

Positive Control (immunohistochemistry): Hippocampus, cerebellum

Species Tested: Human

Related Peptides:

Rattin-24 MAKRGFNCLL LSISEIDLVP KRLE

Rattin-38 MAKRGFNCLL LSISEIDLVP KRLESPNKTR RPYGASIY

Related Antiserum: Anti Rattin-24 Serum Y441

REFERENCES:

- 1)Y. Hashimoto, I. Nishimoto et al., A rescue factor abolishing neuronal cell death by a wide spectrum of familial Alzheimer's disease genes and A β . Proceedings of the National Academy of Sciences 98: 6336-6341, 2001
- 2)H. Tajima, I. Nishimoto et al., Evidence for in vivo production of Humanin peptide, a neuroprotective factor against Alzheimer's disease-related insults. Neuroscience Letters 324: 227-231, 2002

FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



2480-1 Awakura, Fujinomiya-shi,
Shizuoka 418-0011, JAPAN
Phone: +81-544-22-2771
Fax: +81-544-22-2770
URL : www.yanaihara.co.jp
E-mail : ask@yanaihara.co.jp