



## Anti Dioxin Monoclonal Antibody

Cat. No. YII-YM010-EX      Lot No. 5A12B5170531

**Description:** This antibody was raised in a mouse by immunization with porcine thyroglobulin (pTG) conjugate of Dioxin. The product vial contains 1 mg purified IgG<sub>1</sub> of the titled antibody obtained by lyophilizing its 0.01MPBS (pH 7.2, 1 mL) solution without preservative. It can be used for immunoassay or any other immunoreaction with dioxin (2,3,7,8-Tetrachlorodibenzo-*p*-dioxin).

**Immunogen:** 2,3,7,8-Tetrachlorodibenzo-*p*-dioxin-pTG-conjugate      **Host:** Mouse  
**Isotype:** IgG<sub>1</sub> (Protein A purified)

**Product Form:** Lyophilized      **Size:** 1 mg/vial

**Reconstitution:** Reconstitute the product with 1 mL distilled water to make 1 mg/mL stock solution of 0.01M PBS. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN<sub>3</sub> 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antibody solution must be stored at 2°C to 8°C and used within one month. Repeated freezing-thawing should be avoided.

**Suggested Working Dilution Range:** Optimal dilution should be determined by each laboratory for each application.

**Specificity (based on EIA):**

	Dioxins	Crossreactivity
PCDDs	2,3,7,8-Tetrachlorodibenzo- <i>p</i> -dioxin	100.00%
	2,3,7-Trichlorodibenzo- <i>p</i> -dioxin	15.00%
	1,2,3,7,8-Pentachlorodibenzo- <i>p</i> -dioxin	0.60%
	1,2,3,4,6,7,8,9-Octachlorodibenzo- <i>p</i> -dioxin	< 0.10%
PCDFs	2,3,4,7,8-Pentachlorodibenzofuran	14.80%
	1,2,3,4,6,7,8,9-Octachlorodibenzofuran	< 0.01%
Co-PCB	3,3', 4,4', 5-Pentachlorobiphenyl	0.50%
	2,3', 4,4', 5-Pentachlorobiphenyl	< 0.10%

### FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

#### Distributor



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

#### Manufacturer



株式会社  
矢向原研究所

Yanaihara Institute Inc.  
2480-1 Awakura, Fujinomiya-shi,  
Shizuoka 418-0011, JAPAN  
Phone: +81-544-22-2771  
Fax: +81-544-22-2770  
URL : [www.yanaihara.co.jp](http://www.yanaihara.co.jp)  
E-mail : [ask@yanaihara.co.jp](mailto:ask@yanaihara.co.jp)