



Monoclonal Antibody to **Rabbit Smooth Muscle Myosin Heavy Chain(SM1)**

Abbreviation : MCA MH-02

Code No. : 7600

Price : 200 μ l ······ ¥40,000

Clone : 1C10

Quality : BALB/c IgG2b κ

Solution : Ascites (200 μ l)

Character : Working dilution is about 1:3000

Storage & Handling : The solution should be stored below -20°C .

Handling : Repeated freezing and thawing should be avoided.

Note : Research use only.

Reference : 1) Kuro-o, M. et al., J. Biol. Chem. 264, 18272-18275 (1989)

2) Aikawa, M. et al., Circ. Res. 73, 1000-1012 (1993)