



# Monoclonal Antibody to **Nonmuscle Myosin Heavy Chain(SMemb)**

Abbreviation : MCA MH-04

Code No. : 7602

Price : 200  $\mu$ l ······ ¥40,000

Clone : 3H2

Quality : BALB/c IgG2b $\kappa$

Solution : Ascites (200  $\mu$ l)

Character : Working dilution is about 1:3000

Storage & The solution should be stored below  $-20^{\circ}\text{C}$ .

Handling : Repeated freezing and thawing should be avoided.

Note : Research use only.

Reference : 1) Kuro-o, M. et al., J. Biol. Chem. 264, 18272-18275 (1989)

2) Aikawa, M. et al., Circ. Res. 73, 1000-1012 (1993)