

**Monoclonal Antibody to Human Surfactant Protein D (SP-D)**

Code No.: 7608
Quantity: 200 µl (200 µg)
Price: ¥40,000
Clone: 10H11
Quality: BALB/c mouse IgG₁ κ (Purified by protein A chromatography)
Solution: 1mg/ml IgG, 0.1%NaN₃, PBS, pH7.4
Character: This antibody is specific to human SP-D.

Reactivity

to	Human SP-A	<0.2%
	Rat SP-D	<0.2%

Example: This antibody can use for immunohistochemistry. It is possible to react SP-D in tissue section that is deparaffinized and rehydrated on slide.

Storage: The solution should be stored below -20°C. Repeated freezing and thawing should be avoided.

Note: Research use only.

Information**YAMASA CORPORATION****Diagnostics Department**

1-23-8 Nihonbashi Kakigaracho, Chuo-ku, Tokyo 103-0014, Japan

Tel: 81-3-3668-8558 Fax: 81-3-3668-8407

E-mail;dimac@yamasa.com