



Peroxidase Conjugated Monoclonal Antibody to

Human IgG3

Abbreviation: MCA CG-01-HRP
Code No. : 7635
Quantity: 50 μ L (50 μ g)
Clone: 5H9-2a
Quality: BALB/c IgG2b κ -Horseradish Peroxidase
Solution: 1 mg/mL Conjugate, 1% BSA, 0.01% Thimerosal, PBS, pH 7.4

Character:		Reactivity
	to	
	Human IgA	< 0.5%
	Human IgG1	< 0.5%
	Human IgG2	< 0.5%
	Human IgG4	< 0.5%
	Human IgM	< 0.5%

Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.

Note: Research use only



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN

[International/Outside Japan]

TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618

e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com

[国内連絡先]

TEL: 03-5632-9610 / FAX: 03-5632-9619

e-mail: support@cosmobio.co.jp / URL: www.cosmobio.co.jp