



Monoclonal Mouse IgE anti-DNP

Abbreviation: M-IgE
Code No. : 7676
Quantity: 100 μ L (100 μ g)
Clone: 2-9(4)
Quality: Purified from ascites by column chromatographies.
Solution: 1 mg/mL IgE, 0.1% NaN₃, PBS, pH 7.4
Character: This preparation shows single band on SDS-PAGE.

This antibody is specific for DNP.

Storage: The solution should be stored below -20°C.

Repeated freezing and thawing should be avoided.

Note: Research use only

References: 1) Saito, H. et al. J. Immunol. 138, 3927 (1989)

2) Nakamura, H. et.al. Jpn. J. Allergol. 38(12), 1359 (1989)



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN
[International/Outside Japan] [国内連絡先]

TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618

TEL: 03-5632-9610 / FAX: 03-5632-9619

e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com

e-mail: support@cosmobio.co.jp / URL: www.cosmobio.co.jp