



## Human IgE

Abbreviation:	H-IgE
Code No. :	7677
Quantity:	100 $\mu$ L (100 $\mu$ g)
Quality:	Purified by column chromatographies from serum of a IgE myeloma patient. This preparation is negative for HBs antigen, HCV antibody and HIV-1 antibody
Solution:	1 mg/mL IgE, PBS, pH 7.4
Character :	This preparation shows single band on SDS-PAGE. IgG contamination is less than 0.05%.
Storage:	The solution should be stored below -20°C. Repeated freezing and thawing should be avoided.
Note:	Research use only



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN

[International/Outside Japan]

[国内連絡先]

TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618

TEL: 03-5632-9610 / FAX: 03-5632-9619

e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com

e-mail: support@cosmobio.co.jp / URL: www.cosmobio.co.jp