



**Monoclonal Antibody to
Human Cardiac Myosin Light Chain I**

Code No. : 7887
Product name: ANTI-HUMAN CARDIAC MYOSIN LIGHT CHAIN I(520)
Quantity: 200 μ L (200 μ g)
Clone: MLM520
Quality: Balb/c mouse IgG₁ (Purified by protein A chromatography).
Solution: 1mg/mL IgG, 0.1%NaN₃, PBS, pH7.4
Character: This antibody reacts human, porcine, bovine, rat ventricle myosin light chain I.
It is possible to use for western blotting.
Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.
Note: Research use only.