



Peroxidase Conjugated Monoclonal Antibody to Mouse Immunoglobulin Kappa Chain

Code No.: 7962
 Quantity: 200μl(200μg protein)
 Price: ¥30,000
 Clone: 2B6
 Quality: Horseradish peroxidase conjugated Fisher Rat IgG1.
 Solution: 1mg/ml protein, 1% Bovine serum albumin, 0.1% ProClin 300 in PBS, pH7.4
 Storage: The solution should be stored below -20°C.
 Repeated freezing and thawing should be avoided.
 Note: Research use only.
 Character:




Table. Specificities of rat anti-mouse immunoglobulin antibodies

Clone	Code. Purified/HRP Ab	Test clone ^a / isotype								
		1F4 IgG1κ	2B11 IgG1κ	1D1 IgG1λ	5B11 IgG1λ	3H8 IgG2aκ	2G11 IgG2aκ	2F2 IgG2bκ	2A10 IgG2bκ	1D8 IgG3κ
14G11	7965/7966	+	+	+	+	-	-	-	-	-
7D5	7967/7968	-	-	-	-	+	+	-	-	-
2B6	7965/7966	+	+	-	-	+	+	+	+	+
7G12	7963/7964	-	-	+	+	-	-	-	-	-

^a Isotype identified mouse monoclonal antibodies were obtained from Yamasa monoclonal antibody libraries. Anti-mouse immunoglobulin antibodies were applied to mouse monoclonal antibodies-coated ELISA plates. Rat anti-mouse immunoglobulin antibodies that reacted the antigens were probed by horseradish peroxidase-labeled anti-rat IgG and TMB, a peroxidase substrate. +, absorbance > 1.2; - absorbance < 0.15.