



Peroxidase Conjugated Monoclonal Antibody to Mouse Immunoglobulin Lambda Light Chain

Code No.: 7964
Quantity: 200 µl (200 µg protein)
Price: ¥30,000
Clone: 7G12
Quality: Horseradish peroxidase conjugated Fisher Rat IgG2a
Solution: 1mg/ml protein, 1% Bovine serum albumin, 0.1% ProClin 300 in PBS, pH7.4
Storage: The solution should be stored below -20°C.
 Repeated freezing and thawing should be avoided.
Note: Research use only.
Character:




Table. Specificities of rat anti-mouse immunoglobulin antibodies

Clone	Code. Purified/HRP Ab	Test clone ^a / isotype								
		1F4	2B11	1D1	5B11	3H8	2G11	2F2	2A10	1D8
		IgG1κ	IgG1κ	IgG1λ	IgG1λ	IgG2aκ	IgG2aκ	IgG2bκ	IgG2bκ	IgG3κ
14G11	7965/7966	+	+	+	+	-	-	-	-	-
7D5	7967/7968	-	-	-	-	+	+	-	-	-
2B6	7965/7966	+	+	-	-	+	+	+	+	+
7G12	7963/7964	-	-	+	+	-	-	-	-	-

^a Isotype identified mouse monoclonal antibodies were obtained from Yamasa monoclonal antibody libraries. Anti-mouse immunoglobulin antibodies were applied to mouse monoclonal antibodies-coated ELISA plates. Rat anti-mouse immunoglobulin antibodies that reacted the antigens were probed by horseradish peroxidase-labeled anti-rat IgG and TMB, a peroxidase substrate. +, absorbance > 1.2; - absorbance < 0.15.