



Monoclonal Antibody to Rat Surfactant Protein D

Code No. : 80075
Product name: MAB TO RAT SURFACTANT PROTEIN D (5F4)
Quantity: 200 μ L (200 μ g)
Clone: 5F4
Quality: Mouse IgG2a (Purified by protein A chromatography)
Solution: 1 mg/mL IgG, PBS, pH 7.4 (200 μ L)
Character: This antibody reacts with mouse SP-D.
It is possible to use for western blotting and ELISA.
Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.
Note: Research use only.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN
<http://www.cosmobio.com> e-mail : export@cosmobio.co.jp
Phone : +81-3-5632-9617 FAX : +81-3-5632-9618