



Anti Mouse SIRP β Polyclonal Antibody

CatalogNo. GU02-P01AF-A

Lot number	PG-001
Host animal	Rabbit
Immunogen	Recombinant Mouse SIRP β
Purity	Affinity purified
Concentration	0.2 mg/ml
Volume	50 μ l
Formulation	PBS (pH7.4)
Cross Reactivity	Cross reacts with SHPS-1 <ul style="list-style-type: none">• Western Blot: 1/100-1/200. Predicted molecular weight: 50-70 kDa
Application/ Recommended concentration	<ul style="list-style-type: none">• Immunohistochemistry: not tested• Immunoprecipitation: not tested• ELISA: not tested• Flow cytometry: not tested
Storage	Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.
Reference	Akiko Hayashi, Hiroshi Ohnishi, Hideki Okazawa, Seshiru Nakazawa, Hiroshi Ikeda, Sei-ichiro Motegi, Naoko Aoki, Shoji Kimura, Masahiko Mizuki and Takashi Matozaki Positive Regulation of Phagocytosis by SIRP β and Its Signaling Mechanism in Macrophage The Journal of Biological Chemistry Vol.279 No.28, pp29450-29460

For research use only; not for use as a diagnostic.

