



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. GU02-P01AF-A

Anti Mouse SIRP β

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid PBS (pH7.4), Affinity purified
Volume	50 μ l
Concentration	0.2 mg/ml
Specificity	Mouse SIRP β
Antigen	Mouse SIRP β
Isotype	

Application notes WB Other applications are not tested yet.

Recommended use

Recommended dilutions

Western Blot: 1/100-1/200. Predicted molecular weight: 50-70 kDa

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity SHPS-1

Storage Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.

References Hiroshi Ikeda, Sei-ichiro Motegi, Naoko Aoki, Shoji Kimura, Masahiko Mizuki and Takashi Matozaki Positive Regulation of Phagocytosis by SIRP β and Its Signaling Mechanism in Macrophage The Journal of Biological Chemistry Vol.279 No.28, pp29450-29460

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619