



MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. MKM-M05

Anti β 1-receptor, clone 2B9

Product type	Primary antibodies
Host	Mouse
Source	Ascites (Affinity purified)
Form	Liquid PBS (pH 7.4)
Volume	100 μ l
Concentration	1.0 mg/ml
Specificity	β 1 areceptor
Antigen	Synthetic peptide corresponding to the sequence of the second extracellular loop of the human β 1 adrenergic receptor
Clone	2B9
Isotype	IgG2a with κ -type light chains

Application notes WB, ELISA, IP, Immunostaining

Recommended use

Recommended dilutions

Western blotting: 1/100

Immunostaining of human heart tissue and rat primary cultured cardiomyocytes: 1/100

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Human, Rat

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619