



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. MKM-P01

Anti Hyaluronidase

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
Volume	50 µl
Concentration	
Specificity	
Antigen	Synthetic peptide (11 a.a) of human hyaruronidase 2
Isotype	

Application notes

WB, IP

Recommended use

Recommended dilutions

Western blotting: x500

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Human

Storage

Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619