



MONOCLONAL ANTIBODY

Catalog No. PUH-011N-B

Anti Human Cu,Zn-SOD [Super Oxide Dismutase]

Product type	Primary antibodies
Host	Mouse
Source	Culture supernatant
Form	
Volume	500ul
Concentration	1 mg/ml
Specificity	Human Cu,Zn-SOD
Antigen	Human recombinant Cu,Zn-SOD
Clone	5C10
Isotype	IgG1k

Application notes IHC, WB

Recommended use

Recommended dilutions

Western blotting, 1/5,000

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Human

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References


Sugino, N., Shimamura, K.. et al. Changes in activity of superoxide dismutase in the human endometrium throughout the menstrual cycle and in early pregnancy. Human Reproduction, 11, 1073-1078(1996)

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti Human Cu,Zn-SOD [Super Oxide Dismutase]



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619