



MONOCLONAL ANTIBODY

Catalog No. PUH-012N-A

Anti human Cu,Zn-SOD [Super Oxide Dismutase]

Product type	Primary antibodies
Host	Mouse
Source	Culture supernatant
Form	
Volume	100ul
Concentration	1 mg/ml
Specificity	Human Cu,Zn-SOD
Antigen	Human recombinant Cu,Zn-SOD
Clone	5E12
Isotype	IgG2ak

Application notes

WB

Recommended use

Recommended dilutions

Western blotting, 1/5,000

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Human

Storage

Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.


References

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti human Cu,Zn-SOD [Super Oxide Dismutase]



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619