



MONOCLONAL ANTIBODY

Catalog No. PUH-021N-B

# Anti Human Mn-SOD [Super Oxide Dismutase]

<b>Product type</b>	Primary antibodies
<b>Host</b>	Mouse
<b>Source</b>	Culture supernatant
<b>Form</b>	
<b>Volume</b>	500ul
<b>Concentration</b>	1 mg/ml
<b>Specificity</b>	Human Mn-SOD
<b>Antigen</b>	Human Mn-SOD
<b>Clone</b>	2C6
<b>Isotype</b>	IgM

**Application notes** IHC, WB

**Recommended use**

**Recommended dilutions**

Western Blot 1/1,000

Optimal dilutions/concentrations should be determined by the end user.

**Staining Pattern**

**Cross reactivity** Human

**Storage** Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

**References**

Sugino, N., Shimamura, K.. et al. Changes in activity of superoxide dismutase in the human endometrium throughout the menstrual cycle and in early pregnancy. Human Reproduction, 11, 1073-1078(1996)

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**  
Inspiration for Life Science

**Anti Human Mn-SOD [Super Oxide Dismutase]**



←21kDa



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619