



MONOCLONAL ANTIBODY

Catalog No. PUH-051N-A

Anti Hemin

Product type	Primary antibodies
Host	Mouse
Source	Culture supernatant
Form	Liquid
Volume	100ul
Concentration	1 mg/ml
Specificity	Hemin
Antigen	Human gamma-globulin conjugated hemin (naturally prepared)
Clone	1D3
Isotype	IgM

Application notes ELISA
Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage Store below -20°C (below -70°C for prolonged storage).
Aliquot to avoid cycles of freeze/thaw.

References Hifumi, E., Morihara, F., Ishimaru, M., Morikawa, K., Shimizu, K., and Uda, T., Immunochemical Features of Complementarity Determining Region (CDR) Peptide in Anti Hemin Monoclonal Antibody. J. Peptide Sci. 5: 24-31 (1999)

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610
FAX : +81-3-5632-9618 FAX : +81-3-5632-9619