

POLYCLONAL ANTIBODY

Catalog No. TNL-002-SH2

Anti Phospho SHP-1 (Y536)

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
	Filtered (0.22um)
Volume	100 µl
Concentration	
Specificity	
Antigen	GST-mouse SHP-1 (amino acids 516-595, phosphorylated at Y536)
Isotype	IgG

Application notes WB, IHC, IP

Recommended use

Recommended dilutions

Western blotting, 1/1,000 Predicted molecular weight: 68 kDa

Immunohistochemistry, 1/1,000

Immunoprecipitation, 1 µl /1 x 10⁻⁷ cell lysate

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity Rat

Storage Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References 1) Ozawa T, Nakata K, Mizuno K, Yakura H. Negative autoregulation of Src homology region 2-domain-containing phosphatase-1 in rat basophilic leukemia-2H3 cells. Int Immunol. 2007 (9):1049-1061.

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619