

## Anti- PCB #28 (Mouse monoclonal antibody)

---

**Code :** EBT-APM28AS-EX  
**Lot :** S H 0 9  
**Lot size :** 1.0 mL  
**Antigen :** PCB #28-KLH conjugate  
**Isotype :** IgG<sub>1</sub>  
**Purification :** This antibody was isolated from culture supernatant of hybridoma by ProteinG affinity chromatography.  
**Purity :** >80% by SDS-PAGE  
**Concentration :** 0.1 mg /mL  
**Appearance :** Clear liquid  
**Storage :** Short term, 4°C  
Long term, -20°C  
Avoid freeze thaw cycles  
**Preservative :** 0.1% Sodium azide  
**Presentation :** Supplied as a liquid in 0.01M Tris,0.15M NaCl,  
0.1% Sodium azide, pH:7.5  
**Warning :** **For Research Use Only. Not for use in diagnostic procedures.**

---

Distributor



COSMO BIO Co., LTD.  
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618