



Anti-Zebrafish HSP70 : ZHSP67-42 (Mouse monoclonal antibody)

Code : EBT-ZHSP6742-EX
Lot :
Lot size : 1 mL
Isotype : IgG
Antigen : HSP 70 (Zebrafish)
Purity : >80% by SDS-PAGE
Concentration : 0.1 mg /mL
Appearance : Clear liquid
Storage : Short term, 4°C
Long term, -20°C
Avoid freeze thaw cycles
Preservative : 0.1% Sodium azide
Presentation : Supplied as a liquid in 0.01M Tris,0.15M NaCl,
0.1% Sodium azide, pH:7.5

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

e-mail : export@cosmobio.co.jp

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618