



Anti Testosterone-11Alpha-Succinate polyclonal antibody

Product No. FKA-104

The Cross reactivities of anti Testosterone-11Alpha-Succinate (FKA-104) by RI Assay

Antigen	anti Testosterone-11Alpha-Succinate-BSA	
Antibody[FKA-104]	lyophilized Serum	20 μ l / Vial
Dilution of Serum	2000test for RI Assay	250 μ l / tube use
Binding of Total	60-80%	
Testosterone dose of B0/2	150 pg	
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.	

Steroid	Cross Reaction
Testosterone	100 %
5Alpha-Dihydrotestosterone(DHT)	28.50 %
4-Androstenedione	1.55 %
5-Androstene-3Beta,17Beta-diol	0.54 %
Aldosterone	0.07 %
DHA	0.02 %
Progesterone	<0.01 %
Pregnenolone	<0.01 %
17Alpha-Hydroxypregnenolone	<0.01 %
Aldosterone	<0.01 %

For research use only; not for use as a diagnostic.

