



Anti Dihydrotestosterone-11Alpha-Succ polyclonal antibody

Product No. FKA-112E

The Cross reactivities of Anti Dihydrotestosterone-11Alpha-Succ (FKA-112E)
by ELISA with Dihydrotestosterone-11Alpha-Succ-HRP (FKA-111)

Antigen	Dihydrotestosterone-11Alpha-Succ-BSA		
Antibody(FKA-112E)	Purified by Protein A	lyophilized	10µl / Vial
Dilution of IgG	2000test for ELISA	50µl/well use	
Label (FKA-111)	Dihydrotestosterone-11Alpha-Succ-HRP		1ml / Vial
	1,000test for ELISA	50µl/well use	
O.D. (492/600nm) of B ₀	about 0.5		
Testosterone dose of B ₀ /2	about 100 pg		
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.		

Steroid

Cross Reaction

Dihydrotestosterone(DHT)	100 %
Testosterone	58.00 %
4-Androstenedione	10.00 %
Dehydroepiandrosterone	0.02 %
Progesterone	0.03 %
Cortisol	0.01 %
Corticosterone	0.01 %
Estradiol	0.01 %

For research use only; not for use as a diagnostic.

