



Anti Androstenedione-3-CMO polyclonal antibody

Product No. FKA-138

The Cross reactivities of anti Androstenedione-3-CMO (FKA-138) by RI Assay

Antigen	Androstenedione-3-CMO -BSA
Antibody[FKA-138]	lyophilized Serum 10 μ l / Vial
Dilution of Serum	2000test for RI Assay 250 μ l / tube use
Binding of Total	60%
Androstenedione dose of B0/2	about 50 pg/ml
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

Steroid	Cross Reaction
Androstenedion	100 %
5Alpha-Androstanedione	35.00 %
DHA	4.50 %
Androstenetrione	1.50 %
Progesteropne	1.20 %
Testosterone	6.00 %
Cortisol	0.32 %
Corticosterone	0.20 %
Estradiol	0.01 %

For research use only; not for use as a diagnostic.

