



FILCERA™

Asbestos-free Ceramic Fiber Filter Aid

Code No.	: CF101
Appearance	: White granule of short ceramic fiber
Characteristic	: Easily dispersed in water to form excellent filter layer
Purity	: Not less than 99.5 % of alumina and silica
Melting Point	: Not less than 2,000°C
Fiber Diameter	: 2 ~ 4 μm
Fiber Length	: 50 ~ 500 μm
Loss on Ignition	: Not more than 0.05 %
Availability	: 100 g in a box
Use	: Filter aid for quantitative analysis of crude fiber or dietary fiber contained in foods and animal feeds <i>Standard Methods of Analysis in Food Safety Regulation, Japan Food Hygiene Association</i>
References	: <i>Safety Regulation, Japan Food Hygiene Association</i> (1991)

For research use only; not for use as a diagnostic.

