



PRODUCT DATA SHEET

Product: Anti-Amelogenin, Polyclonal

Cat.No: AB 01 (100 μ l/vial)

Species Reactivity:

Bovine ,rat, and mouse. Other species no tested .

Immunogen

Bovine amelogenin .

Storage and Stability:

Store at -20°C or below.

Aliquot to avoid cycles of freeze/thaw

Application and Suggested Dilutions:

Immunohistochemistry:

staining with PAP or ABC methods . Suggested dilution: 100–500 fold .

Western Blot :

Use at a dilution of 1:1,000 to 1:10,000 .

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For in vitro research use only. Not for use in diagnostics or in humans.

Distributor



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



HOKUDO CO.,LTD

4-28,West10,Hachiken9jyo,Nishi-ku,
Sapporo,063-0849,JAPAN
TEL. : +81-11-641-7507
FAX. : +81-11-644-9209
URL : <http://www.hokudo.co.jp>
E-mail : info@hokudo.co.jp