



PRODUCT DATA SHEET

Product: Anti-Gelatinase B

Cat.No: AB36 (100 μ l , 4mg/ml)

Background:

The anti-GelatinaseB polyclonal antibody specifically stains the follicles obtained from the Medaka after ovulation. This antibody may be useful for the study of the influence of environmental hormone-like chemicals on Medaka ovaries.

Host:

Rabbit

Species Reactivity

Medaka, others not tested.

Immunogen

Recombinant Medaka Gelatinase B (MMP9) is a 109-690 amino acid residues and 6xHis.

Format:

100ul of purified rabbit antiserum.

Storage and Stability:

Stable for 2years when stored at less than 4°C.

Aliquot to avoid freeze/thaw cycles.

Application and Suggested Dilutions:

Immunohistochemical Staining :

Incubate frozen sections with antiserum diluted approximately 1:10 for one hour.

Staining may be done with an ACE substrate Kit.

Immunoprecipitation

Western blot

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For in vitro research use only. Not for use in diagnostics or in humans.

Distributor



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer



HOKUDO CO.,LTD

4-28,West10,Hachiken9jyo,Nishi-ku,
Sapporo,063-0849,JAPAN
TEL. : +81-11-641-7507
FAX. : +81-11-644-9209
URL : <http://www.hokudo.co.jp>
E-mail : info@hokudo.co.jp