
Product Data Sheet

Anti HLA class -A,B,C , Mouse Monoclonal antibody

Conjugation : FITC

- Catalog No. : AB-46-F
- Quantity : 100 μ L/vial (1mg/mL)
- Clone : EMR8-5
- Isotype : IgG1
- Buffer : Phosphate buffered saline pH7.4 , 0.05% Sodium Azide
- Immunogen : recombinant HLA-A*2402 extracellular domain
- Specificity : Heavy chains of HLA-A, HLA-B, HLA-C
- Applications : Immunohistochemistry (formalin-fixed tissue)
Western blotting
- Optimal dilutions / concentrations should be determined by the user.
- Storage : Store at 2~8 . Do not freeze.

References :

1. Torigoe T, Asanuma H, Shimozawa K, et al. Immunohistochemical analysis of HLA class I expression in tumor tissues revealed unusually high frequency of down-regulation in breast cancer tissues. submitted
2. Tsukahara T, Kawaguchi S, Torigoe T, et al Prognostic significance of HLA class I expression in osteosarcoma defined by anti-pan HLA class I monoclonal antibody, EMR8-5. *Cancer Sci* 97: 1374-1380, 2006
3. Kitamura H, Torigoe T, Honma I, et al Effect of human leukocyte antigen class I expression of tumor cells on outcome of intravesical instillation of bacillus calmette-guerin immunotherapy for bladder cancer. *Clin Cancer Res* 12: 4641-4644, 2006
4. Komori H, Nakatsura T, Senju S, et al Identification of HLA-A2- or HLA-A24-restricted CTL epitopes possibly useful for glypican-3-specific immunotherapy of hepatocellular carcinoma. *Clin Cancer Res* 12: 2689-2697, 2006
5. Kitamura H, Torigoe T, Honma I, et al Expression and antigenicity of survivin, an inhibitor of apoptosis family member, in bladder cancer: implications for specific immunotherapy. *Urology* 67: 955-959, 2006

For research and in vitro use only. Not for diagnostic or therapeutic work.

HOKUDO CO., LTD

4-28, West 10, Hachiken 9jyo, Nishi-ku, Sapporo, 064-0849, JAPAN

TEL : 81-11-641-7507 FAX : 81-11-644-9209

e-mail : info@hokudo.co.jp URL : <http://www.hokudo.co.jp>