

Product Data Sheet


Anti β 2-Microglobulin Mouse Monoclonal antibody

Catalog No. :	AB47	Lot No. :	GB01B	Exp. :	2008.6
Quantity :	100 μ L/vial (1mg/mL)				
Clone :	EMRB6-12				
Isotype :	IgG1				
Buffer :	Phosphate buffered saline pH7.4 , 0.05% Sodium Azide				
Immunogen :	Recombinant human beta2-microglobulin (B2M)				
Specificity :	human beta2-microglobulin (crossreactive to mouse B2M)				
Applications :	Immunohistochemistry (formalin-fixed tissue , 1:50~1:100) Western blotting				
Storage :	Store at 2~8°C, For long-term storage aliquot and store at -20°C. Avoid freeze / thaw cycles.				
References :	Torigoe T, Asanuma H, Shimosawa K, et al. Immunohistochemical analysis of HLA class I expression in tumor tissues revealed unusually high frequency of down-regulation in breast cancer tissues. J. Translational Med. submitted				

β 2-microglobulin positive



Western Blotting



For research and in vitro use only. Not for diagnostic or therapeutic work.

HOKUDO CO., LTD

4-28, West 10, Hachiken 9jyo, Nishi-ku, Sapporo, 064-0849, JAPAN

TEL : 81-11-641-7507 FAX : 81-11-644-9209

e-mail : info@hokudo.co.jp URL : <http://www.hokudo.co.jp>