



Monoclonal Anti-Human Native α 1-Antitrypsin(AAT)

Technical Data Sheet

Product No.:	IKG-KMSLS3-1
Lot No.:	SA-001
Subclass:	Mouse IgG ₁
Protein(Abs. at 280 nm):	0.5 mg/ml
Quantity:	100 μ g
Clone:	S3-1
Product form:	protein A purified IgG in PBS,0.1%NaN ₃
Storage:	4°C. For long term storage material can be aliquoted and stored at -20°C. Avoid freeze/thaw cycles.
Reactivity:	No cross reactivity with oxidized AAT or elastase-AAT complexes
Use:	EIA,WB

FOR RESEARCH LABORATORY USE ONLY

Distributor

 COSMO BIO CO., LTD.
Inspiration for Life Science
TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME
KOTO-KU, TOKYO 135-0016, JAPAN
TEL : +81-3-5632-9617
FAX : +81-3-5632-9618
URL : <http://www.cosmobio.co.jp/>
e-mail : export@cosmobio.co.jp

Manufacturer

 **Ikagaku Co. Ltd.**
328 Furukawa-cho, Hazukashi, Fushimi-ku,
Kyoto 612-8486 JAPAN
FAX : +81-75-935-0387
E-mail : development@ikagaku.co.jp