



Anti Guinea Pig Fc Receptor Monoclonal Antibody

CATALOG NUMBER:	KEI-FCG27A1
QUANTITY:	200 UG
CLONE:	fcg2-7A1
ANTIGEN:	Guinea Pig Fc Receptor
SOURCE:	Mouse implantation ascites
FORMAT:	Affinity Purified Protein
ISOTYPE:	IgG1
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

