



## Anti Human Calcitonin Monoclonal Antibody

<b>CATALOG NUMBER:</b>	KEI-HCTC08
<b>QUANTITY:</b>	100ug
<b>CLONE:</b>	HCT-C08
<b>ANTIGEN:</b>	Human Calcitonin
<b>SOURCE:</b>	Mouse implantation ascites
<b>FORMAT:</b>	Affinity Purified Protein
<b>ISOTYPE:</b>	IgG1
<b>STORAGE/HANDLING:</b>	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

