



Anti Human Calcitonin Monoclonal Antibody

CATALOG NUMBER:	KEI-HCTC08EX
QUANTITY:	1mg , 10mg
CLONE:	HCT-C08
ANTIGEN:	Human Calcitonin
SOURCE:	Mouse implantation ascites
FORMAT:	Affinity Purified Protein
ISOTYPE:	IgG1
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

