



## Anti Human Beta-Amyloid Monoclonal Antibody

<b>CATALOG NUMBER:</b>	<b>KEI-HMKTB</b>
<b>QUANTITY:</b>	100ug
<b>CLONE:</b>	HMK-TB
<b>ANTIGEN:</b>	Human Beta-Amyloid
<b>SOURCE:</b>	Mouse implantation ascites
<b>FORMAT:</b>	Affinity Purified Protein
<b>ISOTYPE:</b>	IgG1
<b>STORAGE/HANDLING:</b>	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

