



Anti Salmon Calcitonin Monoclonal Antibody

| | |
|--------------------------|--|
| CATALOG NUMBER: | KEI-SCTST01 |
| QUANTITY: | 1mg |
| CLONE: | sCT-ST-01 |
| ANTIGEN: | Salmon Calcitonin |
| SOURCE: | Mouse implantation ascites |
| FORMAT: | Affinity Purified Protein |
| ISOTYPE: | IgG1 |
| STORAGE/HANDLING: | Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles. |

For research use only; not for use as a diagnostic.

