



COSMO BIO Co., LTD.

Inspiration for Life Science

Human Fc epsilon Monoclonal Antibody(CRA1)

CODE NO:	98140
SPECIFICITY:	The monoclonal antibody reacts specifically to alpha-chain of Human Fc epsilon R1
REACTIVITY:	CRA1 reacts with mast cell, basophil and dendritic cell (Langerhans cell).
PROTEIN CONTENT:	Lyophilized powder from ascites containing 200 ug protein. Sterile Filtered. None preservative.
FORM:	This monoclonal antibody was generated by immunizing BALB/C mice and fusing their spleen cells with myeloma cells. Next is derived from ascites fluid and is supplied in PBS. Final are filtered through filters and lyophilized.
MOUSE ISOTYPE:	IgG2b
STORAGE/HANDLING:	After opening aliquot contents and store at -20°C . Store at $2-8^{\circ}\text{C}$ for up to one year. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp