



MOUSE ANTI ALBUMIN MONOCLONAL ANTIBODY

Catalog No.	LNM-KR-001
Clone	A1E8
Concentration	1 mg/ml
Volume	100 UL
Subclass	IgG1k
Immunogen	Human Albumin
Specificity	Human Albumin No cross reactivity with BSA and OVA
Applications	Sandwich ELISA (This antibody can be used as a capture antibody in a human Albumin ELISA in combination with biotinlated, monoclonal antibody A5C9.)
Format	Purified immunoglobulin from BALB/C Mouse Ascites Anmmonium sulfate fractionation Affinity chromatography
Formulation	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
Storage/Handling	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic



抗ヒトアルブミンモノクローナル抗体

品 番	K R - 0 0 1
クローン	A 1 E 8
型	I g G 1, κ
由 来	B A L B / C マウス腹水
精 製 法	硫酸塩析, r - プロテインA
溶 液	P B S (0.09%アジ化ナトリウム含む)
容 量	1 0 0 μ g
保 存 法	4 ~ 8 °C 長期保存の場合は、- 7 0 °C以下にて凍結保存して下さい。 凍結融解は出来るだけ避けて下さい。
備 考	B S A および O V A とは反応しません。 抗ヒトアルブミン抗体 (A1E8) はサンドイッチ ELISA 法の固相に適しています。なお、二次抗体としては抗ヒトアルブミン抗体 (A5C9) が適しています。



コスモ・バイオ株式会社