



COSMO BIO CO., LTD.
Inspiration for Life Science

MOUSE ANTI Ca19.9 MONOCLONAL ANTIBODY

Catalog No.	LNM-KR-013
Clone	2B4
Concentration	1 mg/ml
Volume	100 UL
Subclass	IgG1k
Immunogen	Human CA19.9
Specificity	Purified CA19.9 from tumor cells
Applications	Sandwich ELISA(This one clone is useful for both labeled antibody and capture antibody.)
Format	Purified immunoglobulin from BALB/C Mouse Ascites Anmmonium sulfate fractionation Affinity chromatography
Formulation	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
Storage/Handling	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp

抗ヒトCa19.9モノクローナル抗体

品 番	K R - 0 1 3
クローン	2 B 4
型	I g G 1, κ
由 来	B A L B / C マウス腹水
精 製 法	硫酸塩析, r - プロテインA
溶 液	P B S (0.09%アジ化ナトリウム含む)
容 量	1 0 0 μ g
保 存 法	4 ~ 8 °C 長期保存の場合は、- 7 0 °C以下にて凍結保存して下さい。 凍結融解は出来るだけ避けて下さい。
備 考	抗ヒトCa19.9抗体(2B4)はサンドイッチELISA法に適しています。



コスモ・バイオ株式会社