



MOUSE ANTI ELASTASE 1 MONOCLONAL ANTIBODY

Catalog No.	LNM-KR-017
Clone	207E10
Concentration	1 mg/ml
Volume	100 UL
Subclass	IgG1k
Immunogen	Human Elastase 1
Specificity	Human Elastase 1
Applications	Sandwich ELISA (This antibody can be used as a capture antibody in a human Elastase 1 ELISA in combination with biotinlated, monoclonal antibody 203F4.)
Format	Purified immunoglobulin from BALB/C Mouse Ascites Anmmonium sulfate fractionation Affinity chromatography
Formulation	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
Storage/Handling	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic

