



MOUSE ANTI ELASTASE 1 MONOCLONAL ANTIBODY

Catalog No.	LNM-KR-018
Clone	203F4
Concentration	1 mg/ml
Volume	100 UL
Subclass	IgG2bk
Immunogen	Human Elastase 1
Specificity	Human Elastase 1
Applications	Sandwich ELISA (This antibody can be used as a detection antibody in a human Elastase 1 ELISA combination with the monoclonal capture antibody 207E10.)
Format	Purified immunoglobulin from BALB/C Mouse Ascites Anmmonium sulfate fractionation Affinity chromatography
Formulation	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
Storage/Handling	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic



抗ヒトエラスターゼ 1 モノクローナル抗体

品 番	K R - 0 1 8
クローン	2 0 3 F 4
型	I g G 2 b , κ
由 来	B A L B / C マウス腹水
精 製 法	硫酸塩析, r - プロテイン A
溶 液	P B S (0.09% アジ化ナトリウム含む)
容 量	1 0 0 μ g
保 存 法	4 ~ 8 °C 長期保存の場合は、- 7 0 °C 以下にて凍結保存して下さい。 凍結融解は出来るだけ避けて下さい。
備 考	抗ヒトエラスターゼ 1 抗体 (203F4) はサンドイッチ ELISA 法の二次抗体に適しています。なお、固相としては抗ヒトエラスターゼ 1 抗体 (207E10) が適しています。



コスモ・バイオ株式会社